



Vascular access in the ICU

Understanding the role of PICCs in the critical care setting



Meeting critical patient needs

Selecting an appropriate vascular access device (VAD) in the intensive care unit (ICU) is complex, meriting careful examination of infusion needs and other considerations.

Infusion needs

In the ICU, critically ill patients often require central venous access with a range of performance requirements, including:

- Rapid intravenous infusion of fluids for rehydration purposes¹
- Intravenous infusion of drugs and solutions that are vesicants, irritants, and/or are hyperosmotic¹
- Administration of parenteral nutrition¹
- Need for repeated blood sampling¹
- Hemodynamic monitoring¹



Other considerations

When selecting central venous access, consider these factors:

- Patient diagnosis and comorbidities¹
- Potential complications during placement and removal¹
- Recommended patient positioning¹
- Projected time in the ICU¹
- Anticipation of post-ICU therapy¹



Catheters in the ICU: PICCs, FICCs, and CICC

In the ICU, there are a variety of catheter options to meet critical patient needs. Carefully consider which may be best for each situation.

		Power injectable PICCs	Centrally Inserted Central Catheters (CICCs)	Upper Groin Femorally Inserted Central Catheters (FICCs)
	Terminate in a central vein	✓ ²	✓ ³	✓ ⁴
	Deliver incompatible infusates via multi-lumen catheter	✓ ⁵	✓ ³	✓ ⁶
	Allow for blood draw	✓ ²	✓ ³	✓ ⁷
	Pressure-injectable for use with CT scans	✓ ²	✓ ⁸	✓ ⁸
	Compatible with CVP monitoring	✓ ²	✓ ³	✓ ⁷
	Recommended for long-term use	✓ ²	✗ ³	✗ ⁷
	Recommended use beyond 48 hours	✓ ²	✓ ⁸	✗ ⁷
	Can be placed and removed while patient is slightly elevated	✓ ⁹	✗ ³	✓ ¹⁰
	No risk of carotid, subclavian, or femoral artery puncture during insertion	✓ ¹¹	✗ ¹¹	✗ ¹¹
	No risk of iatrogenic pneumothorax	✓ ²	✗ ¹²	✓ ¹³
	No risk of iatrogenic hemothorax	✓ ²	✗ ¹⁴	✓ ⁶
	Removal does not require supine or Trendelenburg	✓ ²	✗ ³	✓ ¹⁵
	Accommodates flow rates for rapid resuscitation	✗	✓ ¹⁶	✓ ⁶

Advantages of placing PICCs in the ICU

- Central access for high-acuity complex cases such as collars, burns, halos²
- No risk of iatrogenic pneumothorax and iatrogenic hemothorax²
- Reduced risk of bleeding during placement compared to other central lines¹⁷
- Microbial load demonstrated to be lower in peripheral insertion sites in comparison to neck and chest¹⁸
- More options for patient positioning during insertion and removal^{2,9}
- Reduced need for confirmatory X-ray with tip confirmation technology compared to blind placement¹⁹
- Can be used for short- and long-term vascular access²

Clinical advantages of placing PICCs in the ICU




PICCs can provide central access for higher-acuity patient scenarios

When a patient enters the ICU, consider the patient's condition, potential for long-term treatment, and overall health and comfort.

Examples of conditions favoring PICCs:²⁴

- Anatomical (helmets, halo brace, collars)
- Coagulopathies (clotting disorders, anticoagulant drugs)
- Severe pulmonary conditions (respiratory failure, included COVID-19 pneumonia)
- Tracheostomy
- Elevated risk of pleuropulmonary complications
- Burns or infected skin in the chest/neck area
- Obesity



PICCs carry no risk of iatrogenic pneumothorax and hemothorax²

When selecting an insertion site for central venous access, carefully consider the risks of insertion location. There are several life-threatening CICC insertion-related risks that are not associated with PICCs.^{14,20}


Iatrogenic Pneumothorax

An iatrogenic pneumothorax is an unexpected hospital-acquired condition where there is the presence of air or gas in the cavity between the lungs and the chest wall, causing collapse of the lung.²¹

Iatrogenic Hemothorax

An iatrogenic hemothorax is an unexpected hospital-acquired condition where blood collects in the space between the chest wall and the lung (the pleural cavity).²²

Rates for image-guided central venous access		
Insertion Site Location	Iatrogenic Pneumothorax	Iatrogenic Hemothorax
Subclavian	2%-3% ²³	1% ²³
Jugular	1% ²³	1% ²³
FICC	0% ²³	0% ²³
PICC	0% ²³	0% ²³



The risk of pneumothorax and hemothorax when using CICCs is low; however, it is still present.²³

Iatrogenic pneumothorax and hemothorax is associated with subclavian and jugular catheter insertions.^{12,23} For patients with lung conditions—such as pulmonary conditions, lung cancer, and/or pneumonia—it's important to minimize the risk of additional complications (like iatrogenic pneumothorax and/or hemothorax) as compromised respiratory systems are at play before insertion risks are introduced. A PICC may provide the least amount of risk of adverse complications of iatrogenic pneumothorax and hemothorax during insertion.^{12,24}




PICCs are associated with reduced risk of major arterial puncture compared to other central lines

It is worth considering the arteries near various catheters insertion sites and the patient's risk of bleeding (coagulopathy).

Consider this:

Vessels commonly utilized for insertion sites of CICCs, PICCs, and FICCs	Vessels in close proximity
Internal jugular vein ²⁷	Carotid artery ²⁷
Subclavian vein ²⁷	Subclavian artery ²⁷
Femoral vein ²⁷	Femoral artery ²⁸
Brachial vein ²⁹	Brachial artery ²⁹
Basilic vein ²⁸	Brachial artery ²⁸



Consider the Guidelines on Coagulopathy

According to the Society of Interventional Radiology (SIR) and endorsed by the Cardiovascular and Interventional Radiological Society of Europe (CIRSE), PICCs have a “lower risk of bleeding, being more easily detectable and controllable.”¹⁷ This is a Category 1 recommendation.

Researchers have additionally noted that “arterial access in the periphery is easily managed with digital compression-promoting coagulation.”³⁰

While ultrasound is a widely available tool that can aid in proper catheter placement, researchers have also noted that the use of ultrasound does not eliminate the risk of internal jugular arterial puncture.^{27,31}



Chest and neck have a higher colonization rate than peripheral access sites³²

Did you know?

According to the ESCAPE³ study, **microbial load** at the **catheter site impacts risk** of increased colonization during

- Insertion of vascular access devices (VAD)¹⁸
- Long-term dwell of the device¹⁸


Also that **site selection**

- May have an important role in preventing a potential harmful infection¹⁸

¹⁸Evaluation of skin colonization and placement of vascular access device exit sites—a study to determine if presence of varying levels of microorganisms on the skin had an effect on vascular access device (VAD) site selection.¹⁸



The relationship of microbial load on skin at vascular access insertion site



Higher microbial load

Femoral and internal jugular access sites are associated with increased risk of central line-associated bloodstream infection (CLABSI)³² because they have a higher microbial load on the skin during and after device placement

Lowest microbial load

was found in the upper arm in the largest percent of population of patients in the ICU¹⁸

Know the risks: tracheostomies and infection

When a patient has a tracheostomy, they may be at higher risk for central line infection.³³ It's recommended to stay far away from the jugular site when placing a central line; therefore, a PICC may be a more appropriate choice in this scenario.

Site selection—The American Society of Anesthesiologists recommends selecting a catheter insertion site that is not contaminated or potentially contaminated (i.e. adjacent to a tracheostomy).³³

Risk—The Society for Healthcare Epidemiology of America states that patients who have a tracheostomy concurrent with a jugular site catheter placement may be at higher risk for infection.³⁴





PICCs offer more options for patient positioning during insertion and removal

Patient positioning for placement and removal of central lines			
	CICCs	FICCs	PICCs
Trendelenburg	✓ ³⁵	✗	✗
Prone	✗	✗	✓ ³⁶
Supine	✓ [*]	✓ ⁶	✓ ⁹
Semi-Fowlers	✗	✓ ⁷	✓ ⁹

*Trendelenburg recommended for CICC placement and removal by the Joint Commission; supine if necessary

Consider this:

For placement

According to recent studies,³⁷ it's been shown that PICCs can be inserted while the patient is sitting up.

CICCs can only be inserted when the patient is in Trendelenburg, or supine if necessary,³⁵ according to the Joint Commission.

** Gruppo Accessi Venosi Centrali Alungo Termine

For removal

PICCs No published cases of air embolism associated with PICC removal have been reported¹⁵

CICCs & FICCs 30 mins

INS guidelines recommend that patient should be flat for 30 mins to reduce air embolism and bleeding risk after removal of CICC and FICC.¹⁵



PICCs can be placed at the bedside with navigation and tip confirmation³⁸

The Infusion Nurses Society recommends to “position the tip of a CVAD in the lower third of the superior vena cava (SVC) at or near the cavoatrial junction (CAJ) for adults and children”.¹⁵ With both **bedside navigation and tip confirmation technology, PICCs can help with first-attempt success rates and correct tip positioning**¹⁹ — so you can confidently place the catheter and help optimize patient comfort.

Procedural placement technology support	BD power injectable PICCs ²	CICC ³⁹	FICC ⁴⁰
Able to be placed with navigation systems	✓	✓ [*]	✓ [*]
Usable with tip confirmation systems	✓	✓ [*]	✓ [*]
Able to be inserted with ultrasound	✓	✓	✓

*Dependent on manufacturer of device, length of catheter, and complimentary confirmation system.

It's acceptable to use a PICC line as soon as you place it when a p-wave is present, reducing the need for X-ray confirmation.⁴¹





PICCs can be used for short-term and long-term treatment

When a patient enters the ICU, consider the patient's condition, potential for long-term treatment, and overall health and comfort.



FICC

FICC use is limited to emergent placement by these authoritative organizations:

- Infusion Nursing Society¹⁵
- The Joint Commission⁴²
- Asia Pacific Society of Infection Control⁴³

If a FICC or other emergent line is placed under suboptimal aseptic technique, it should be replaced within 24 to 48 hours.^{15,43}




CICC

Manufacturers may limit their CICCs to short-term use. CICCs may require insertion of multiple catheters throughout the patient's hospital stay resulting in more sticks, corresponding with high risk of central line associated bloodstream infection.³



PICCs are indicated for short and long-term dwell and can be used throughout the continuum of care.^{19,26}



Putting patient assessment first

When a patient is admitted to the ICU, placing a PICC may **reduce insertion-related complications** (compared to other central vascular access devices), which may lead to mortality and morbidity.²⁰ PICC placement is **not dependent on supine or Trendelenburg positioning**.^{9,36} PICCs offer an alternative vascular access solution for ICU patients to central lines inserted through the neck, chest, or groin and can **follow the patient throughout the continuum of care**.

For assistance on finding vascular access management solution for your patients and hospital, **speak to your BD representative**.



References

- 1 M. Scoppettuolo G. *The GAVeCeLT Manual of Picc and Midline: Indication, Insertion, Management*. 1st ed. Edra S.p.A. 2020.
- 2 C. R. Bard, Inc. Polyurethane PICC with Microintroducer. Bard Access Systems Inc; 2016.
- 3 ARROWgard Blue Plus® antimicrobial multi-lumen central venous catheterization product with sharps safety features; K-45703-108A (3-02). Arrow International; 2022.
- 4 Hanlon G, on behalf of the ANZICS Safety and Quality Committee. *ANZICS Central Line Insertion and Maintenance Guideline*. 2012.
- 5 McPherson DJ. Topics in progressing care: Peripherally inserted central catheters: What you should know. *Nursing 2008 Critical Care*. 2008;3(3):10-15.
- 6 Pal N. Femoral central venous access periprocedural care. *Medscape*. Published July 5, 2022. Accessed July 12, 2022. <https://emedicine.medscape.com/article/80279-periprocedure>
- 7 Castro D, Martin Lee LM, Bhutta BS. Femoral vein central venous access. In: *StatPearls*. Treasure Island, FL: StatPearls Publishing; September 12, 2022.
- 8 ARROW® Pressure Injectable CVC; 2011-0212v2. Arrow International; 2022.
- 9 Mitsuda S, Tokumine J, Matsuda R, Yorozu T, Asao T. PICC insertion in the sitting position for a patient with congestive heart failure: A case report. *Medicine (Baltimore)*. 2019;98(6):e14413.
- 10 Annetta MG, Marche B, Dolcetti L, et al. Ultrasound-guided cannulation of the superficial femoral vein for central venous access. *J Vasc Access*. 2022;23(4):598-605. doi:10.1177/11297298211003745
- 11 Leib AD, England BS, Kiel J. Central Line. In: *StatPearls*. Treasure Island (FL): StatPearls Publishing; Updated July 25, 2022.
- 12 Tsotsolis N, Tsigogianni K, Kioumis I, et al. Pneumothorax as a complication of central venous catheter insertion. *Ann Transl Med*. 2015;3(3):40. doi:10.3978/j.issn.2305-5839.2015.02.11
- 13 Shieh L, Go M, Gessner D, Chen JH, Hopkins J, Maggio P. Improving and sustaining a reduction in iatrogenic pneumothorax through a multifaceted quality-improvement approach. *J Hosp Med*. 2015;10(9):599-607. doi:10.1002/jhm.2388
- 14 Hashimoto Y, Fukuta T, Maruyama J, Omura H, Tanaka T. Experience of peripherally inserted central venous catheter in patients with hematologic diseases. *Intern Med*. 2017;56(4):389-393. doi:10.2169/INTERNALMEDICINE.56.7625
- 15 Gorski LA, Hadaway L, Hagle ME, et al. Infusion therapy standards of practice. *J Infus Nurs*. 2021;44(suppl 1):S1-S224. doi:10.1097/NAN.0000000000000396
- 16 Roe III EJ. Central venous access via infraclavicular (subclavian/subclavicular) approach to subclavian vein. *Medscape*. Published October 17, 2021. Accessed July 12, 2022. <https://emedicine.Medscape.com/article/80336>
- 17 Patel IJ, Davidson JC, Nikolic B, et al. Consensus guidelines for periprocedural management of coagulation status and hemostasis risk in percutaneous image-guided interventions. *J Vasc Interv Radiol*. 2012;23(6):727-736. doi:10.1016/j.jvir.2012.02.012
- 18 Moureau N, Marsh N, Zhang L, et al. Evaluation of skin colonisation and placement of vascular access device exit sites (ESCAPE study). *J Infect Prev*. 2019;20(1):51-59. doi:10.1177/1757177418805836
- 19 Kwon S, Son SM, Lee SH, et al. Outcomes of bedside peripherally inserted central catheter placement: a retrospective study at a single institution. *Acute Crit Care*. 2020;35(1):31-37. doi:10.4266/acc.2019.00731
- 20 Pittiruti M, Brutti A, Celentano D, et al. Clinical experience with power-injectable PICCs in intensive care patients. *Crit Care*. 2012;16(1):R21. doi:10.1186/cc11181
- 21 Ojeda Rodriguez JA, Hipskind JE. Iatrogenic Pneumothorax. In: *StatPearls*. Treasure Island (FL): StatPearls Publishing; November 13, 2021.
- 22 Borke J. Hemothorax: Medlineplus medical encyclopedia. *MedlinePlus*. Published July 1, 2020. Accessed June 23, 2022. <https://medlineplus.gov/ency/article/000126.htm>
- 23 Dariushnia SR, Wallace MJ, Siddiqi NH, et al. Quality improvement guidelines for central venous access. *J Vasc Interv Radiol*. 2010;21(7):976-981. doi:10.1016/j.jvir.2010.03.006
- 24 Pinelli F, Little A, Kokotis K, Alsbrooks K, Pittiruti M. Assessment of the MAGIC recommendations in context of evolving evidence based on use of PICC in the ICU. *J VAS Access*. 2021;11297298211048019. doi:10.1177/11297298211048019
- 25 Camp-Sorrell D, Matey L. *Access Device Standards of Practice for Oncology Nursing*. Pittsburgh: Oncology Nursing Society. 2017.
- 26 Ryu DY, Lee SB, Kim GW, Kim JH. A peripherally inserted central catheter is a safe and reliable alternative to short-term central venous catheter for the treatment of trauma patients. *J Trauma Inj*. 2019;32(3):150-156.
- 27 Troianos CA, Hartman GS, Glas KE, et al. Guidelines for performing ultrasound guided vascular cannulation: recommendations of the American Society of Echocardiography and the Society of Cardiovascular Anesthesiologists. *J Am Soc Echocardiogr*. 2011;24(12):1291-1318.
- 28 Safety Committee of Japanese Society of Anesthesiologists. Practical guide for safe central venous catheterization and management 2017. *J Anesth*. 2020;34(2):167-186.
- 29 O'Grady NP, Alexander M, Burns LA, et al. Summary of recommendations: guidelines for the prevention of intravascular catheter-related Infections. *Clin Infect Dis*. 2011;52(9):1087-1099. doi:10.1093/cid/cir138
- 30 Hill S. Avoiding complications during insertion. In: *Vessel Health and Preservation: The Right Approach for Vascular Access*. Moureau N (ed). Springer Cham. 2019. <https://doi.org/10.1007/978-3-030-03149-7>
- 31 Ruetzler K, Hauffe T, Piegeler T, Hatem A, Spahn DR. An accidental arterial puncture and anticoagulation after internal jugular vein catheterization resulting in massive hematoma and airway compromise. *Int J Anesthetic Anesthesiol*. 2018. doi:10.23937/2377-4630/1410069
- 32 Mollee P, Jones M, Stackelroth J, et al. Catheter-associated bloodstream infection incidence and risk factors in adults with cancer: a prospective cohort study. *J Hosp Infect*. 2011;78(1):26-30. doi:10.1016/j.jhin.2011.01.018
- 33 Practice guidelines for central venous access 2020: An updated report by the American Society of Anesthesiologists Task Force on Central Venous Access. *Anesthesiology*. 2020;132(1):8-43. doi:10.1097/ALN.0000000000002864
- 34 Marschall J, Mermel L, Fakih M, et al. Strategies to prevent central line-associated bloodstream infections in acute care hospitals: 2014 update. *Infect Control Hosp Epidemiol*. 2014;35(7):753-771. doi:10.1086/676533
- 35 The Joint Commission. *Preventing Central Line-Associated Bloodstream Infections: Useful Tools, An International Perspective*. Nov 20, 2013. Accessed July 12, 2022. <http://www.jointcommission.org/CLABSIToolkit>
- 36 Patrona-Aurand R, Lovett KA, Kochur K. Innovative positioning for peripherally inserted central catheter insertion on a prone patient. *J Assoc Vasc Access*. 2016;21(4):212-216. doi:10.1016/j.java.2016.06.004
- 37 Narita A, Takehara Y, Maruchi Y, et al. Usefulness of peripherally inserted central catheter port system (PICC-PORT) implantation in the sitting position: a new technique for cases unsuitable for conventional implantation. *Jpn J Radiol*. 2022;8. doi:10.1007/s11604-022-01317-7
- 38 Becton, Dickinson and Company. (2021). Economic and Clinical Value of Bedside PICC Placement.
- 39 Brescia F, Pittiruti M, Ostroff M, Spencer TR, Dawson RB. The SIC protocol: a seven-step strategy to minimize complications potentially related to the insertion of centrally inserted central catheters. *J Vasc Access*. 2021;11297298211036002. doi:10.1177/11297298211036002
- 40 Brescia F, Pittiruti M, Ostroff M. Rapid femoral vein assessment (RaFeVA): A systematic protocol for ultrasound evaluation of the veins of the lower limb, so to optimize the insertion of femorally inserted central catheters. *J Vasc Access*. 2020;00(-):1-10. doi:10.1177/1129729820965063
- 41 Gao Y, Liu Y, Zhang H, Fang F, Song L. The safety and accuracy of ECG-guided PICC tip position verification applied in patients with atrial fibrillation. *Ther Clin Risk Manag*. 2018;14:1075-1081.
- 42 The Joint Commission. Oakbrook Terrace, IL: Joint Commission Resources; May 2012.
- 43 Ling ML, Apisarnthanarak A, Jaggi N, et al. APSIC guide for prevention of Central Line Associated Bloodstream Infections (CLABSI). *Antimicrob Resist Infect Control*. 2016;5:16. doi:10.1186/s13756-016-0116-5

Becton, Dickinson and Company
Salt Lake City, UT 84116, USA
800.545.0890

bd.com

