




# Animal Health Catalogue

Anaesthesia Systems



## BD™ Quincke spinal needle

- Key/slot arrangement of stylet and cannula hubs facilitates aligned needle and stylet bevel orientation
- Translucent window hub (featuring contact clarity) that helps visualization of cerebrospinal fluid
- Fitted stylet designed to reduce tissue coring
- Color-coded
- Sterile, single use



Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box/ Shipping Case	Reference
27 G	3 1/2"	90	Grey	25	405259
26 G	3 1/2"	90	Brown	25	405258
25 G	4 11/16"	119	Blue	10	405234
25 G	3 1/2"	90	Orange	25	405257
25 G	3"	75	Orange	25	405246
25 G	2"	51	Blue	25	405245
25 G	1"	25	Blue	25	405243
23 G	3 1/2"	90	Turquoise	25	405240
22 G	7"	178	Black	10	405149
22 G	5"	127	Black	10	405148
22 G	3 1/2"	90	Black	25	405256
22 G	3"	75	Black	25	405255
22 G	2 1/2"	63	Black	25	405244
22 G	1 1/2"	38	Black	25	405254
20 G	6"	152	Yellow	10	405211
20 G	3 1/2"	90	Yellow	25	405253
20 G	3"	75	Yellow	25	405252
20 G	1 1/2"	38	Yellow	25	405251
19 G	3 1/2"	90	Cream	25	405250
19 G	3"	75	Cream	25	405249
18 G	3 1/2"	90	Pink	25	405248
18 G	3"	75	Pink	25	405247

## BD™ Quincke spinal needle

With Introducer 22 G 1 1/4"

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
27 G	3 1/2"	90	Grey	25	405069

With Introducer 20 G 1 1/4"


Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
26 G	3 1/2"	90	Brown	25	405065
25 G	3 1/2"	90	Orange	25	405084

BD™ Quincke Spinal Needle introducers

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
20 G	1 1/4"	32	Yellow	25	405260
22 G	1 1/4"	32	Black	25	405261

## BD Whitacre™ pencil point spinal needle

- Manufactured with proprietary needle tipping, side-port and low-glare hub technologies
- Tip designed to spread the dural fibres and help reduce the occurrence of post dural puncture headache
- The thin wall design is intended for improved visualization of cerebrospinal fluid
- Pencil point design intended to provide better “tip strength” to minimize bending or breakage
- Translucent hub for ease of cerebrospinal fluid visualization
- Color-coded
- Sterile, single use



Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
27 G	3 1/2"	90	Grey	25	402051
27 G	4 11/16"	119	Grey	10	409443
25 G	3 1/2"	90	Blue	25	402050
25 G	4 11/16"	119	Blue	10	409442
24 G	3 1/2"	90	Purple	25	405104
22 G	3 1/2"	90	Black	25	401995

## BD Whitacre™ pencil point spinal needle with introducer sets

### With Integral Introducer 22 G 1 1/4"

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
27 G	3 1/2"	90	Grey	25	405075
27 G	4 1/16"	103	Grey	25	405113

### With Integral Introducer 20 G 1 1/4"

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
25 G	3 1/2"	90	Blue	25	405076
25 G	4 1/16"	103	Blue	25	405112

## BD™ Tuohy and BD™ Tuohy Weiss epidural needle

- Epidural needles with locking base and 'modified Tuohy' bevel intended to minimize accidental dural puncture
- Needle graduation marks every 10 mm, paediatric epidural needle graduation marks every 5 mm
- Needle equipped with moulded plastic hub
- Color-coded
- Sterile, single use



## Tuohy Epidural Needles

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
17 G	3 1/2"	90	Grey	10	405092
18 G	3 1/2"	90	Pink	10	405018
20 G	3 1/2"	90	Yellow	10	405028

## Tuohy Weiss Epidural Needles winged

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
16 G	3 1/2"	90	White	10	405194
17 G	3 1/2"	90	Grey	10	405191
17 G	5"	127	Grey	10	405193
18 G	3 1/2"	90	Pink	10	405188
18 G	5"	127	Pink	10	405190
20 G	3 1/2"	90	Yellow	10	405186



## BD Epilor™ Loss-of-Resistance (L.O.R) Syringes

Description	Shipping Case	Reference
7 ml L.O.R. Plastic Syringe with BD Luer-Lok™ Tip	10	405291
7 ml L.O.R. Plastic Syringe Luer Slip Tip	10	405292

## Kits for epidural anaesthesia procedures with BD™ Weiss needle

A selection of kits for epidural anaesthesia procedures with BD™ Weiss needles.

- Modified Tuohy needle point designed for greater penetration characteristics
- Shorter point length than standard epidural needles designed to minimize accidental dural puncture
- BD Epilor™ Loss-of-Resistance (L.O.R.) plastic syringe designed for better detection of epidural space
- Sterile, single use



### Perisafe™ Plus! – Epidural Kit

Description	Shipping Case	Reference
1 Epidural Tuohy Weiss needle (with wings) 17 G x 80 mm 1 Closed end Epidural Catheter 19G x 36" (80 mm) 1 loss-of-resistance syringe (L.O.R.) 7 ml 1 Flat filter 0.22 µ 1 Connector	12	401506M
1 Epidural Tuohy Weiss needle (with wings) 18 G x 80 mm 1 Closed end Epidural Catheter 20G x 36" (914 mm) 1 loss-of-resistance syringe (L.O.R.) 7 ml 1 Flat filter 0.22 µ 1 Connector	12	401508M

### Mini-Perisafe™ Plus! – Epidural Kit

Description	Shipping Case	Reference
1 Epidural Tuohy Weiss needle (with wings) 18 G x 80 mm 1 Epidural Catheter 20 G x 36" (914 mm) 1 Flat filter 0.22 µ 1 Connector 1 BD Luer-Lok™ cap	12	401511M

### BD Perisafe™ – Epidural Mini-Kit

Description	Shipping Case	Reference
1 Epidural Tuohy Weiss needle 18 G x 90 mm 1 Closed end Epidural Catheter 20 G x 36" (914 mm) 1 loss-of-resistance syringe (L.O.R.) 7 ml 1 Flat filter 0.22 µ 1 Connector 1 BD Luer-Lok™ cap	10	400273

BD Switzerland Sàrl, Terre Bonne Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland  
Tel: +41 21 556 30 00

**bd.com**

© 2019 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company. CF04131