



The leader in hypodermic products

Hypodermic product catalog

A complete line of medication administration and safety-engineered products for the healthcare industry



Table of contents

BD SafetyGlide™ Syringes and Needles (sterile, single use)	2
BD Eclipse™ Syringes and Needles	3
BD Integra™ Syringes and Needles	3
BD PrecisionGlide™ Needles (sterile, single use)	4
BD® Blunt Fill/Filter Needles	4
BD® Syringes (sterile, single use)	5
BD® Syringe/Needle Combination (sterile, single use)	6
BD® Allergy Syringes	7
BD® Tuberculin Syringes	7
BD® Oral and Enteral Syringes	8

BD SafetyGlide™ Syringes and Needles (sterile, single use)

BD SafetyGlide™ Hypodermic Needle

50 per box, 500 per case

Cat no.	Description
305921	27 G x 5/8 in.
305916	25 G x 1 in.
305901	25 G x 5/8 in.
304387	23 G x 1 1/2 in.
305902	23 G x 1 in.
305917	21 G x 1 1/2 in.
305915	21 G x 1 in.
305918	18 G x 1 1/2 in.



BD SafetyGlide™ Needle with Syringe

50 per box, 400 per case

Cat no.	Description
1 mL BD Luer Slip Syringe with Detachable SafetyGlide™ Needle	
305927	27 G x 5/8 in.
305903	25 G x 5/8 in.
3 mL BD Luer Slip Syringe with Detachable SafetyGlide™ Needle	
305924	25 G x 1 in.
305904	25 G x 5/8 in.
303326	23 G x 1 1/2 in.
305905	23 G x 1 in.
305909	21 G x 1 1/2 in.



Scan this code to watch a video on how to use the BD SafetyGlide™ Needle:



1 mL BD SafetyGlide™ Tuberculin Syringe with Permanently Attached Needle

100 per box, 400 per case

Cat no.	Description
305945	27 G x 1/2 in. regular bevel
303327	27 G x 3/8 in. intradermal bevel



1 mL BD SafetyGlide™ Allergy Syringe Tray with Permanently Attached Needle

25 per tray, 1,000 per case

Cat no.	Description
305950	27 G x 1/2 in. regular bevel
303330	27 G x 3/8 in. intradermal bevel



BD Eclipse™ Syringes and Needles

BD Eclipse™ Needle

100 per box, 1,200 per case

Cat no.	Description
305757	30 G x 1/2 in.
305758	27 G x 1/2 in.
305838	25 G x 1 1/2 in.
305837	25 G x 1 in.
305829	25 G x 5/8 in.
303304	23 G x 1 1/2 in.
305826	23 G x 1 in.
305765	21 G x 1 1/2 in.
305764	21 G x 1 in.
305712	18 G x 1 1/2 in.

Scan this code to watch a video on how to use the BD Eclipse™ Needle:



BD Eclipse™ Needle and Syringe

50 per box, 300 per case

Cat no.	Description
1 mL BD Luer-Lok™ Syringe with Detachable Eclipse™ Needle	
305778	30 G x 1/2 in.
305789	27 G x 1/2 in.
305780	25 G x 5/8 in.
3 mL BD Luer-Lok™ Syringe with Detachable Eclipse™ Needle	
305787	25 G x 1 in.
305781	25 G x 5/8 in.
303306	23 G x 1 1/2 in.
305782	23 G x 1 in.
305784	21 G x 1 1/2 in.
305779	21 G x 1 in.
5 mL BD Luer-Lok™ Syringe with Detachable Eclipse™ Needle	
303307	23 G x 1 1/2 in.
10 mL BD Luer-Lok™ Syringe with Detachable Eclipse™ Needle	
303321	21 G x 1 1/2 in.



BD Integra™ Syringes and Needles

3 mL BD Integra™ Syringe with Detachable Needle

100 per box, 400 per case

Cat no.	Description
305270	25 G x 1 in.
305269	25 G x 5/8 in.
305271	23 G x 1 in.
305274	21 G x 1 1/2 in.
305273	21 G x 1 in.

Scan this code to watch a video on how to use the BD Integra™ Syringe:




BD PrecisionGlide™ Needles (sterile, single use)

BD® Conventional Needles

100 per box, 1,000 per case

Cat no.	Description
305128	30 G x 1 in.
305106	30 G x 1/2 in.
301629	27 G x 1 1/2 in.
305109	27 G x 1/2 in.
305110	26 G x 3/8 in. intradermal bevel
305127	25 G x 1 1/2 in.
305125	25 G x 1 in.
305122	25 G x 5/8 in.
305194	23 G x 1 1/2 in.
305145	23 G x 1 in.
305129	21 G x 2 in.
305167	21 G x 1 1/2 in.
305165	21 G x 1 in.
305196	18 G x 1 1/2 in.
305195	18 G x 1 in.
305198	16 G x 1 1/2 in.
305197	16 G x 1 in.

Types of syringe tips



BD Luer-Lok™ Tip sizes:

1 mL
3 mL
5 mL
10 mL
20 mL
30 mL
50 mL

Slip tip sizes:

1 mL
3 mL
5 mL
10 mL
20 mL
30 mL
50 mL

Eccentric tip sizes:

20 mL
50 mL


Catheter tip sizes:

50 mL

Types of needle walls and bevels

The thin wall needle is thinner than a regular wall needle. The inner diameter is bigger without any increase in the outside diameter. For example, a 21 G thin wall needle has the same outside diameter as a 21 G regular wall, but has the same inner diameter as a 20 G needle. Thin wall allows the same flow volume as a regular wall needle one gauge larger. Thin wall is available in 16, 18, 21, 23 and 25 G. Regular wall needles are available in 26, 27 and 30 G.

Regular wall




Thin wall

Regular bevel



Intradermal bevel



BD® Blunt Fill/Filter Needles

100 per box, 1,000 per case

Cat no.	Description
305180	18 G x 1 1/2 in. BD® Blunt Fill Needle
305181	18 G x 1 in. BD® Blunt Fill Needle
305211	18 G x 1 1/2 in. BD® Blunt Filter Needle (5 micron filter)
305201	18 G x 1 1/2 in. BD Nokor™ Filter
305213	16 G x 1 in. BD Nokor™ Vented
305214	18 G x 1 in. BD Nokor™ Vented
305215	18 G x 1 1/2 in. BD Nokor™ Admix
305216	16 G x 1 in. BD Nokor™ Admix

Scan this code to watch a video on how to use the BD® Blunt Fill and BD® Blunt Filter Needles:



BD® Conventional Needle and Syringe color chart:



BD[®] Syringes (sterile, single use)

1 mL BD Luer-Lok[™] Syringe only

100 per box, 800 per case (1/100 mL)

Cat no.	Description
309628	BD Luer-Lok [™] Tip



3 mL BD[®] Syringe only

200 per box, 800 per case (1/10 mL)


Cat no.	Description
309657	BD Luer-Lok [™] Tip
309656	Slip tip



5 mL BD[®] Syringe only

125 per box, 500 per case (1/5 mL)


Cat no.	Description
309646	BD Luer-Lok [™] Tip
309647	Slip tip



10 mL BD[®] Syringe only

100 per box, 400 per case (1/5 mL)

Cat no.	Description
302995	BD Luer-Lok [™] Tip
303134	Slip tip



25 per box, 100 per case (1/5 mL)

Cat no.	Description
309695	BD Luer-Lok [™] Tip Control Syringe


20 mL BD[®] Syringe only

48 per box, 192 per case (1 mL)

Cat no.	Description
302830	BD Luer-Lok [™] Tip
302831	Slip tip

120 per box, 480 per case (1 mL)

Cat no.	Description
300613	Eccentric tip



30 mL BD[®] Syringe only

56 per box, 224 per case (1 mL)

Cat no.	Description
302832	BD Luer-Lok [™] Tip
302833	Slip tip




50 mL BD[®] Syringe only

40 per box, 160 per case (1 mL)

Cat no.	Description
309653	BD Luer-Lok [™] Tip
309654	Slip tip

60 per box, 240 per case (1 mL)


Cat no.	Description
300866	Eccentric tip



50 mL BD[®] Catheter Tip Syringe only (with tip shield)

40 per box, 160 per case

Cat no.	Description
309620	Catheter tip



BD[®] Pharmacy Convenience Pak



Cat no.	Description	Tray contents
309701	1 mL BD Luer-Lok [™] Slip Tip Syringe	25 syringes per sterile tray, 12 trays per case of 300
309702	3 mL BD Luer-Lok [™] Tip Syringe	25 syringes per sterile tray, 12 trays per case of 300
309703	5 mL BD Luer-Lok [™] Tip Syringe	25 syringes per sterile tray, 12 trays per case of 300
309605	10 mL BD Luer-Lok [™] Tip Syringe	20 syringes per sterile tray, 12 trays per case of 240
305617	20 mL BD Luer-Lok [™] Tip Syringe	10 syringes per sterile tray, 12 trays per case of 120
305618	30 mL BD Luer-Lok [™] Tip Syringe	10 syringes per sterile tray, 12 trays per case of 120
309680	50 mL BD Luer-Lok [™] Tip Syringe	20 syringes per sterile tray, 6 trays per case of 120



BD[®] Syringe/Needle Combination (sterile, single use)

3 mL BD[®] Syringe/Needle Combination BD Luer-Lok[™] Tip

100 per box, 800 per case (1/10 mL)

Cat no.	Description
309582	25 G x 1 1/2 in.
309581	25 G x 1 in.
309570	25 G x 5/8 in.
309589	23 G x 1 1/2 in.
309571	23 G x 1 in.
309577	21 G x 1 1/2 in.
309575	21 G x 1 in.
309580	18 G x 1 1/2 in.
305060	18 G x 1 1/2 in. BD [®] Blunt Fill Needle



5 mL BD[®] Syringe/Needle Combination BD Luer-Lok[™] Tip

100 per box, 400 per case (1/5 mL)

Cat no.	Description
309633	21 G x 1 1/2 in.
309632	21 G x 1 in.
305062	18 G x 1 1/2 in. BD [®] Blunt Fill Needle



10 mL BD[®] Syringe/Needle Combination BD Luer-Lok[™] Tip

100 per box, 400 per case (1/5 mL)

Cat no.	Description
309643	21 G x 1 1/2 in.
309642	21 G x 1 in.
305064	18 G x 1 1/2 in. BD [®] Blunt Fill Needle



BD Luer-Lok[™] Sterile Syringe Tip Cap



200 per box, 2,000 per case

Cat no.	Description
305819	BD [®] Syringe Tip Cap

BD Luer-Lok[™] Sterile Syringe Tip Cap Tray

10 per tray, 200 per box, 2,000 per case

Cat no.	Description
308341	BD [®] Syringe Tip Cap Tray



BD[®] Allergy Syringes

1 mL BD SafetyGlide™ Allergy Syringe Tray with Permanently Attached Needle

25 per tray, 1,000 per case



Cat no.	Description
305950	27 G x 1/2 in. regular bevel
303330	27 G x 3/8 in. intradermal bevel

1/2 mL BD[®] Allergy Tray with Permanently Attached Needle

25 per tray, 1,000 per case



Cat no.	Description
305535	27 G x 1/2 in. regular bevel
305536	27 G x 3/8 in. intradermal bevel

1 mL BD[®] Allergy Tray with Permanently Attached Needle

25 per tray, 1,000 per case



Cat no.	Description
305540	27 G x 1/2 in. regular bevel
305541	27 G x 3/8 in. intradermal bevel
305542	27 G x 3/8 in. regular bevel

BD[®] Tuberculin Syringes

1/2 mL BD[®] Tuberculin Syringe with Permanently Attached Needle

100 per box, 500 per case (1/100 mL)


Cat no.	Description
305620	27 G x 1/2 in. regular bevel
300860	27 G x 3/8 in. intradermal bevel



1 mL BD[®] Tuberculin Syringe with Detachable Needle, slip tip

100 per box, 800 per case (1/100 mL)


Cat no.	Description
309623	27 G x 1/2 in. regular bevel
303344	27 G x 3/8 in. intradermal bevel
309626	25 G x 5/8 in. regular bevel
309624	21 G x 1 in. regular bevel



1 mL BD[®] Tuberculin Syringe only

200 per box, 1,600 per case (1/100 mL)

Cat no.	Description
309659	Slip tip



BD® Oral and Enteral Syringes

BD® Clear Oral Syringe

100 per box, 500 per case

Cat no.	Description
305217	1 mL
305220	3 mL
305218	5 mL
305219	10 mL

BD® Amber Oral Syringe

100 per box, 500 per case

Cat no.	Description
305207	1 mL
305210	3 mL
305208	5 mL
305209	10 mL



BD® Enteral Syringe with UniVia™ Connection

Cat no.	Description	Packaging
305851	1 mL	200 per box, 1,600 per case
305853	3 mL	200 per box, 800 per case
305855	5 mL	125 per box, 500 per case
305857	10 mL	100 per box, 400 per case
302837	20 mL	48 per box, 192 per case
302836	30 mL	56 per box, 224 per case
305863	50 mL	40 per box, 160 per case



BD® Enteral Syringe Tip Cap

100 per box, 1,000 per case

Cat no.	Description
302435	BD® Enteral Syringe Tip Cap



BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

[bd.com](https://www.bd.com)

BD and the BD Logo, Autoguard, Eclipse, Integra, Luer-Lok, PrecisionGlide and SafetyGlide are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2023 BD. All rights reserved. BD-75999 (03/23)

