

# BD Pyxis™ ES Refrigerator

Secure, single-line item storage of high-risk and high-value refrigerated medication

## What fits in a standard/half-height bin

Medication type*	Quantity*
Vials 1 mL	175
Vials 2 mL	100
Vials 4 mL	75
Vials 5 mL	70
Vials 10 mL	24
Syringes 1 mL	50
Syringes 2 mL	30
Syringes 4 mL	20
Syringes 4 mL plus needles in separate pkg.	10
Syringes 10 mL	12
Ampules 1 mL	40
Suppositories	75
Eye drops 15 mL bottles	30
Rabies vaccine kits	4

\*Medication packaging may vary by manufacturer. Represented quantity is an approximation.



Standard/half-height bin internal dimensions (*w x h x d*) 4.85 in. x 1.52 in. x 12.57 in.



## What fits in a deep/full-height bin

Medication type*	Quantity*
Syringes 10 mL	40
PCA syringes 30 mL	16
PCA syringes 40/50 mL	10
IV bags 50 mL (with overwrap packaging)	6
IV bags 50 mL (without overwrap packaging)	20
IV bags 100 mL (with packaging)	4
IV bags 100 mL (without packaging)	10
IV bags 250 mL	8
IV bags 500 mL	4
IV bags 1 liter	2
Bottles 250/300 mL	4
Bottles 500 mL	2

\*Medication packaging may vary by manufacturer. Represented quantity is an approximation.



Deep/full-height bin internal dimensions (w x h x d) 4.85 in. x 4.52 in. x 12.57 in.

## Common types of refrigerated medications

- Insulin
- Vaccines (e.g., hepatitis, influenza, rabies)
- Antibiotics (e.g., penicillin)
- IV solutions
- Suppositories
- Neuromuscular blockers and paralytics (e.g., Rocuronium/Zemuron™, succinylcholine chloride/Quelicin®)
- Lorazepam/Ativan®
- Compounded IV solutions
- Oncology drugs (e.g., octreotide/Sandostatin®, TBO-filgrastim/Granix™)
- Famotidine/Pepcid®
- Alteplase/Cathflo® Activase®
- Darbepoetin alfa/Aranesp®
- Oral antibiotic suspensions (e.g., amoxicillin clavulanate/Augmentin®, cephalexin/Keflex®)

BD, San Diego, CA, 92130, U.S.

[bd.com](http://bd.com)

BD, the BD Logo and Pyxis are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2019 BD. All rights reserved. 191-1691 (0319/PDF) BD-9725.

